

STOOL COLLECTION INSTRUCTIONS

If the request is for . . .

Stool Culture
Enteric Pathogens
Salmonella/Shigella
C&S
Enterohemorrhagic E. coli
(eg: *E. coli* 0157)



Transport Media
Used by CLS:

C&S Media
(Orange Cap)

Store at Room Temp.

If the request is for . . .

Stool for Parasites
Ova & Parasite
Giardia
Amoeba
Intestinal Parasites
Worms
Giardia Antigen
Cryptosporidium Antigen



Transport Media
Used by CLS:

Para-Pak
10% Buffered
Neutral
Formalin
(Pink cap)

If the request is for . . .

C. difficile Toxin/Culture
Fecal Leukocytes
Fecal WBC
Rotavirus



Transport Media
Needed - Refrigerate.

Clean Transport
Container

Store at 2°-8°C (refrigerate)



CONSOLIDATED
LABORATORY SERVICES
St. Vincent's HealthCare

For all your laboratory needs
308-5600

REV: 8/14